
Otto Bock®

QUALITY FOR LIFE

4R148

(D)	Harmony® e-pulse	3
(GB)	Harmony® e-pulse	10
(F)	Harmony® e-pulse	17
(I)	Harmony® e-pulse	24
(P)	Harmony® e-pulse	31
(NL)	Harmony® e-pulse	38
(S)	Harmony® e-pulse	45
(DK)	Harmony® e-pulse	52
(N)	Harmony® e-pulse	59
(TR)	Harmony® e-pulse	66

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Please read these instructions for use carefully and pay special attention to the safety considerations!

1 Components

The Harmony® e-pulse is delivered with the following components

- | | |
|---------------------------|----------------|
| 1. Harmony® e-pulse | 4R148 |
| 2. Strap Attachment Plate | 4X226 |
| 3. Screws (2x) | 501S138=3.0×16 |
| 4. Torx Key | 710H11=10IP |
| 5. Filter | 4X221 |
| 6. Charger | 757L16-2 |
| 7. Plug Adapter | 4X219 |
| 8. Socket Connector | 2R117=0 |

2 Safety Considerations

Caution!

If the Harmony® e-pulse or any of the accompanying parts show any signs of damage, the correct operation of the device cannot be guaranteed and the device must not be used!

Caution!

Avoid all explosive/combustible materials while using the Harmony® e-pulse!

Notice!

Do not use the Harmony® e-pulse without the filter that is included in the delivery. Dust particles can damage the pump!

Notice!

Only charge the Harmony® e-pulse when the prosthesis is not being worn!

Notice!


Protect the Harmony® e-pulse from moisture! Do not flush the pump with water! Do not immerse in water!

Notice!

Protect the Harmony® e-pulse from any mechanical damages!

Notice!

The Harmony® e-pulse must be charged and used within the recommended environmental conditions. Extreme heat and cold can damage device functions!

3 Description

3.1 Application

The Harmony® e-pulse is to be used for prosthetic fittings. It is a component of the Harmony® System.

The Harmony® e-pulse does not have a system height, and thus is particularly well-suited for transfemoral or transtibial patients with a long residual limb.

3.2 Indication

Indications according to the Otto Bock Mobility System MOBIS®:

Intended for patients of MOBIS Grade 1 – 4; Weight classes 1—4.

3.3 Effect

The Harmony® e-pulse is an electronic pump that can create an elevated vacuum in a prosthetic socket, thereby improving the connection between the amputee and the prosthesis. The system consists of a gel liner, an airwick, a rigid socket (which is connected to the pump), a sealing sleeve, and the Harmony® e-pulse pump. The pump creates an elevated vacuum in the space of the airwick between the rigid socket and liner when the Harmony® e-pulse is activated.

System Advantages:

- Reduces daily volume differences
- Improved adhesion/suspension
- Improved proprioception
- Reduces the forces that are applied to the residual limb

4 Application

Caution!

The socket has to be fabricated following the Harmony® socket fabrication technique. The use of other techniques may cause health risks for the patient!

We recommend attaching the Harmony® e-pulse to the socket using the 4X226 Strap Attachment Plate (Fig. 2). This plate can be attached with the provided attachment screws 501S138=3.0×16 by using the Torx Key 710H11=10IP. The screw threads in the housing are made of plastic: Do not overtighten the screws!

Once the plate is in place, the Harmony® e-pulse can be attached with a strap to the socket, the cosmesis or the pylon. The shorter the distance between pump and Socket Connector the quicker the desired vacuum can be achieved.

The Harmony® e-pulse can also be mounted upside down. This does not influence pump function.

To connect the Harmony® e-pulse with the Socket Connector 2R117 shorten the Silicon in-take tube (Fig. 1: IN) to the required length and use the 4X221 Filter as a connector. If a cosmesis is used, the expulsion tube (Fig. 1: OUT) can be routed through the foam to the outside.

Note: Shortening the expulsion tube could increase the noise level of the pump.

5 Functions

5.1 Basic Functions

The following Harmony® e-pulse functions can be operated by the patients themselves:

Function	Action
On	Push + Button, short
Off	Hold - Button, for at least 2 seconds. Note: While holding the - Button the pump valve is open and air can stream into the socket. The longer the button is being pressed the more air can stream back.
Adjust Vacuum Level	Push either the + or - Button. The + Button increases the vacuum level (Sound: short, high), while the - Button reduces the vacuum level each time it is pressed (Sound: short, medium). The LEDs glow according to the chosen level. Notes: • When in sleep mode, the LED lights will not be lit. To wake device out of sleep mode, press any button once. The LED lights will activate and show current status. Once the e-pulse has woken up, pressing any button will begin adjustments. • Any time the Harmony® e-pulse has been reset it will start up in Level 1 when first activated. If the system has not been reset, it will always start up in the previously used level. • The LED next to the plug socket indicates level 1 (minimum), the one next to the Otto Bock logo indicates level 4 (maximum). • When the minimum or maximum adjustment is reached a short, deep tone will be heard. • If the pump cannot reach the requested vacuum level within a certain time frame the Harmony® e-pulse will make a short warning tone and the 2nd LED light will begin flashing. If this occurs it is very likely that the system has a leak.
Message Mode	To activate: Hold + Button for 5 seconds. To de-activate: Push any button or simply move the system—a built-in movement sensor will recognize the movement and switch the device into normal mode. Note: The Message Mode alternately increases and decreases the vacuum level.
Sound on/off	Push + and - Buttons simultaneous. Note: Tones accompany the different functions/actions. See the table below for the meaning of the various tones.
Charging	Remove the plug socket dust cover and insert the plug. The charging process will be accompanied by an increasing LED signal. To indicate that the charging process is complete, all four LED lights will glow simultaneously. Advice: • Only charge the Harmony® e-pulse when the prosthesis is not being worn. • Depending on the vacuum level, the socket volume or activity level, the battery can power the system for 2-3 days per charge. Nevertheless daily charging is recommended. • The battery should be charged once per day (best during the night) to ensure that the system will operate all day long. The battery capacity and battery life will not be affected by daily charging. But nevertheless the capacity depends on the battery age and may decrease over the life time of the Harmony® e-pulse.
Reset	If the device does not appear to be working properly, reset the system by holding + and - Buttons simultaneously for approx. 5 seconds. When the device starts up again it will automatically default to level 1.

The following chart gives an overview about the LED and warning tones for each function:

Symbol Overview

Harmony® e-pulse

- Off
- Glowing (Status or Adjustment)
- ⚡ Flashing (Error or Warning)

Function	Action	Visual Feedback (LEDs)	Description	Audible Feedback	Tone/Length
On	Push Button	According to last used vacuum level			
Off	Hold Button		Fading out		Falling/Long
Increase Vacuum	Push Button		Increasing		High/Short
Decrease Vacuum	Push Button		Decreasing		Medium/Short
Massage Mode on	Hold Button (5 Sec.)		Changing		
Massage Mode off	Any Button or Movement				
Sound on/off	Push & Button				
Reset	Hold & Button (5 Sec)				High/Long

Signals	Visual Feedback (LEDs)	Description	Audible Feedback	Tone/Length
Maximum/Minimum vacuum setting reached				Low/Short
Battery is charging		Increasing, beginning with LED next to plug		
Battery charge completed		All 4 LEDs glowing		

5.2 Low-Power-Mode

The Low-Power-Mode is an optional feature, which is not activated on delivery. This feature should only be adjusted by the prosthetist who does the patient fitting.

The Harmony® e-pulse has an internal movement sensor. It allows the Harmony® e-pulse to automatically recognize resting phases and switch itself to the energy saving Low Power Mode. In the Low-Power-Mode the vacuum level will slowly be reduced to a lower level. Any incipient movement or any Button push will cause the system to automatically reactivate to the user level.

To activate the Low Power Mode follow the following steps:

Low-Power-Mode	Action
Activate Low Power Mode	1. Bring the device in to charging state by applying the AC charger. 2. Hold  Button for approx. 10 seconds. 3. Adjust to the lowest requested level by pushing either the  or  Button. 4. To save the adjustments hold  Button for approx. 10 seconds until the device executes a reset. The Low Power Mode is now activated. Unplug the charger. Notes: • When entering the adjustment menu for the first time (as delivered) all LED lights will glow when the mode is deactivated. • When re-entering the menu the LED lights will show the current adjustment status. If all LED lights glow it means that the mode is deactivated. • If there has been no adjustment for 10 seconds the Harmony® e-pulse automatically exits this menu and returns to the original setting.
Deactivate Low Power Mode	1. Bring the device in to charging state by applying the AC charger. 2. Hold  Button for approx. 10 seconds. 3. Turn on all LED lights by pushing the  Button. 4. To save the adjustments hold  Button for approx. 10 seconds. The Low Power Mode is now deactivated. Unplug the charger.

6 Maintenance

The Harmony® e-pulse should be examined at least once per year to check for any wear or loss of functionality.

You may need to replace the external filter 4X221 when

- a) it shows strong colour changes and/or
- b) the Harmony® e-pulse seems to pump normally but does not create an elevated vacuum between socket and liner. To check the filter: Remove the Socket Connector tube from the filter and activate the pump. If the pump stops despite the open end: Change the filter.

Tip!

If the pump operates constantly or too often while wearing the prosthesis and a leak has been excluded, try to „air flush“ the system. Take the prosthesis off and let the pump operate a few times until the time limit function stops it. This way dust particles may be transported out of the system. Never use water to flush the pump!

7 Warning and Error Signals

Signal	Action
Low battery	The battery will be empty soon. The amount of time remaining can vary depending on socket volume, activity, outside temperature etc. and the set vacuum level. The Harmony® e-pulse should be charged as soon as possible following the instructions under 5.1 "Charging". If there is no opportunity to charge the device the sound can be switched off as described under 5.1 "Sounds on/off". The warning then only reoccurs as a short vibration.
Empty battery	The battery is empty and the Harmony® e-pulse switches off. Recharge the Harmony® e-pulse according to step 5.1"Charging".
Battery Error	This error may occur if the system is getting too hot (e.g. because of the charging process or direct exposure to the sun, etc.). Remove the device from the charger or from the source of heat, switch it off and let it cool down. If the error still occurs when the cooled-down device is switched on, the Harmony® e-pulse needs to be returned to Otto Bock for service.
Leak in the system	The Harmony® e-pulse has a timer to help determine how long the pump needs to achieve a certain elevated vacuum level. If the preset time limit is exceeded the pump will give the user a warning to indicate that the system has a probable leak. The first thing to check should always be the proximal seal between liner and socket. If this does not solve the problem check the tube and filter connections.

Symbol Overview
Harmony® e-pulse

○ Off

● Glowing (Status or Adjustment)

⚡ Flashing (Error or Warning)

Warnings and Errors	Visual Feedback (LEDs)	Description	Audible Feedback	Tone/Length
Battery low	○ ⚡ ⚡ ○	Changing		High-Low series/Single beeps
				If sound off: 3 x Vibrations
Battery empty - Device turning off	● ○ ○ ○			Falling/Long
Battery error - Battery overheating	⚡ ○ ○ ●			High-Low siren
Leak in the system	○ ⚡ ○ ○			High-Low siren

8 Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the Instructions for Use. The device must only be operated with tested modular components in accordance with the Otto Bock Mobility System. The manufacturer is not liable for damage caused by component combinations that were not authorized by the manufacturer.

9 CE conformity

This product meets the requirements of the 93/42/EEG guidelines for medical products. This product has been classified as a class I product according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by Otto Bock with sole responsibility according to appendix VII of the guidelines.

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Otto Bock HealthCare GmbH